



Please move icons and buttons as marked

“MP1, MP2,...MP5” correspond to:

- MP1 = LIQUID_PUMP1,
- MP2 to LIQUID PUMP2
- MP3 to LIQUID PUMP3
- MP4 to GAS_PUMP1
- MP5 to GAS_PUMP2
- and they work like microbio app’s buttons: if Manual/Auto == Auto, then the buttons show in GREEN/RED for ON of OFF but buttons cannot be clicked. If == Manual, they also show GREEN/RED for ON/OFF but now they are also clickable buttons that change state of the pumps if clicked. When changing from AUTO to manual as “Are you sure you want to turn off auto mode”?